# Trace40C<sup>®</sup>

PRODUCT SHEET

## LANGUAGES

**p. 3**      **ENGLISH**

p. 9      ITALIANO

## THE PROBLEM

Ovarian cancer (OC) is one of the most lethal cancers in women, with 300,000 cases diagnosed worldwide in 2018 and 180,000 deaths.

The lack of accurate screening and diagnostic tools and the rapid spread of disease are the main causes of lethality.

Biological markers, such as Cancer Antigen-125 (CA-125) have been shown highly sensitive but poorly specific for screening, and qualitative transvaginal ultrasound (US) imaging are not enough informative to detect OC at the early stage in the general population.

The failure of a screening program leads to often useless surgeries for benign masses at final histology and to the lack of centralization of oncological cases with worsening of the patient's prognosis. An adequate and reproducible preoperative diagnosis is therefore of paramount importance.


## OUR SOLUTION: TRACE4OC®

The adnexal ovarian masses appear to have characteristics at ultrasonography (US) that are helpful in the diagnosis to expert operators. However, classification guidelines developed so far are helpful in high and low risk group discrimination, while showing inconclusive for lesions assigned with an intermediate risk.

TRACE4OC is our proprietary licensed software that uses machine learning to deeply analyze US images of adnexal ovarian masses of women at risks, in combination with the CA-125 level and with the current woman's menopausal state to provide the individual level of risk (very low risk or medium-high risk) of being affected by ovarian cancer.

## DESCRIPTION

**TRACE4OC provides the woman's level of risk (very low or medium-high risk) of being affected by ovarian cancer, by a semi-automatic reading of the woman's ovarian mass, obtained from a transvaginal US image frame, without contrast medium, in combination with CA-125 measure and the menopausal state of the woman.**

## **INDICATION OF USE of TRACE4OC®**

TRACE4OC is indicated for woman at risk of ovarian cancer on the recommendation of gynecology specialists.

Physicians with expertise with pelvic US can use TRACE4OC to support the reporting on US-detected adnexal masses investigations in these subjects. Specialists in gynecology can use TRACE4OC as an aid to diagnosis.

TRACE4OC software provides the woman level of risk (very low risk or medium-high risk) of being affected with ovarian cancer at the date of the TRACE4OC-processed US investigation of an adnexal mass, combined with the serum level test for CA-125, performed no earlier than and no later than one month from the US study, and with the current woman's menopausal state (premenopause or postmenopause).

However, it should be noted that DeepTrace Technologies S.R.L. considers TRACE4OC as a support to gynecology specialists in their diagnosis, and as a support to physicians with expertise with pelvic US in their reporting on adnexal masses. These medical professionals keep the final decision-making responsibility.

## HOW TRACE4OC® PREDICTS THE RISK OF OVARIAN CANCER

TRACE4OC is a stand-alone software based on US image analytics in the ovarian region, data mining and machine learning classifiers. It is based on the identification of US image features specifically associated to abnormal morphometric and echogenic pattern within and near an ovarian adnexal cancer mass caused by the ovarian cancer disease, in combination with the serum level test for CA-125, corrected for the menopausal state (pre-menopause, post-menopause).

With TRACE4OC, ovarian mass features of a woman are extracted from her US pelvis image (obtained by transvaginal probe without contrast medium), then automatically analyzed by algorithms of support vector machines and random forests, two families of machine learning classifiers that, once properly trained, allow automatically distinguishing between two classes of women of interest, in this case:

**-women with medium-high risk of ovarian cancer vs**

**-women with very low risk of ovarian cancer.**

TRACE4OC uses the features of the morphometric and echogenic pattern of the adnexal ovarian masses of the classes of women described above, masses of which the histological diagnosis is known, to train its learning machines to distinguish the images of the two classes. TRACE4OC training is supervised by the histological diagnosis of the masses (malignant disease, benign disease) obtained following definitive surgery.

When the US image of the ovarian mass of a woman at risk for ovarian cancer is provided as input to TRACE4OC, in combination with the CA-125 measurement and the woman's menopausal status, the software uses its learning machines to classify that ovarian mass into one of two classes (medium-high risk of ovarian cancer vs very low risk of ovarian cancer).

## PERFORMANCES

Performance of TRACE4OC was tested on a set of US transvaginal images obtained from patients of the AROMA/MULTIAROMA study<sup>1</sup> with suspicious ovarian masses, known CA-125 serum level and menopausal state.

Overall, in the AROMA study, TRACE4OC was used to predict the risk of malignancy of the ovarian masses based only on the US appearance, performing with a weighted mean of 84% accuracy, 79% sensitivity, and 86% specificity (Chiappa V, J. Ultrasound, 2020).

In the MULTIAROMA study, TRACE4OC was used to predict the risk of malignancy of the ovarian masses based on the US appearance, on the acoustic shadow presence, on the woman's level of serum CA-125 and menopausal state, showing a weighted mean of 89% accuracy, 99% sensitivity, 77% specificity (Chiappa V, Eur rad Exp, 2021).

*Example of US transvaginal image of an ovarian mass for a patient at risk of ovarian cancer, with highlighted the segmented region of interest of the ovarian mass.*

**ACCURACY 89%\***  
**SENSITIVITY 99%**  
**SPECIFICITY 77%**

\* performance obtained using US in combination with CA-125 and menopausal state. The performance obtained using only US are: ACCURACY 84%, SENSITIVITY 79%, SPECIFICITY 86%

## SCALABILITY

The present version of TRACE4OC has been tailored using US transvaginal ovarian images of subjects at risks of ovarian cancer, but the training set can be customized for other US modalities including US transabdominal, as well as biological and clinical biomarkers. The software can be customized depending on users' needs and available data.

[1] Protocol MULTIAROMA Istituto Nazionale dei Tumori (Milan, Italy) N157/20 at the date of 30 April 2020.

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